
# APXMED

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Ningbo Greatcare Trading Co., Ltd., Unit 93, Building 12, No. 818 Qiming Road, Yinzhou, 315105 Ningbo, Zhejiang, China (Китай).

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## Contents

---

- 23 Abbreviation list
- 23 Symbols
- 24 Introduction
- 24 Product designation
- 24 General description
- 26 Safety information
- 27 Direction for use
- 33 Care and maintenance
- 34 Storage and transportation
- 34 Disposal information
- 34 Warranty

## Abbreviation list

mm Hg	millimeter of mercury
°C	degree Celsius
%	percent
©	copyright

## Symbols

See user manual

Keep dry

Temperature limit

Keep away from sunlight

Contains latex

Should be disposed separately from the normal waste stream

Date of manufactured

Catalogue number

Batch code

## Introduction

---

Apexmed International BV wishes to thank you for purchasing one of its ranges of products and you are invited to read this user manual thoroughly. This user manual contains all of the information necessary to use the device properly. Please adhere closely to the notices contained within and keep the user manual in a place where it is easy to find and is accessible to the operator. The contents of this manual can be changed, without warning or additional obligations, therefore, the manufacturer can make changes and improvements. It is prohibited to copy or translate any of this user manual without the written permission of Apexmed International B.V.

---

**Attention!** Before using the appliance read the user manual.

---

## Product designation

---

Aneroid blood pressure device is intended for determination of blood pressure by Korotkoff's method. The device can be used in clinics as well as by individuals.

## General description

---

The device consists of diaphragm manometer, compression cuff with connection tubes and pneumatic inflation bulb. AT-12, AT-13 и AT-22 models are also packed with stethoscopes.

## The work principle

Blood pressure measuring principle is based on listening to Korotkoffs sounds with stethoscope in decompression mode on distal area of brachial artery that is preliminary pinched by compression cuff overpressure, and determination of the pressure level on diaphragm manometer scale.

## Variants

- AT-10 — cuff with metal fixing ring, without stethoscope, with diaphragm pressure gauge 50 mm;
- AT-11 — cuff without metal fixing ring, without stethoscope, with diaphragm pressure gauge 50 mm;

- AT-12 — cuff with metal fixing ring and stethoscope built in a cuff, with manometer 50 mm;
- AT-13 — cuff with metal fixing ring, with stethoscope, with manometer 50 mm;
- AT-20 — combined PALM type with cuff, metal fixing ring and diaphragm pressure gauge 65 mm, without stethoscope;
- AT-21 — combined PALM type with cuff, without metal fixing ring, with diaphragm pressure gauge 65 mm, without stethoscope;
- AT-22 — combined PALM type, with cuff, metal fixing ring and diaphragm pressure gauge 65 mm, stethoscope built in a cuff
- AT-30 — combined PALM type with cuff, with metal fixing ring, with diaphragm pressure gauge 50 mm, without stethoscope;
- AT-31 — combined PALM type with cuff, without metal fixing ring, with diaphragm pressure gauge 50 mm, without stethoscope;
- AT-41 — Desktop model with cuff and with large face, without metal fixing ring, without a stethoscope.

### Technical characteristics

Measuring range, mm of mercury (cuff pressure)	0–300
Admissible error, mm of mercury (overall measuring range)	±3
Air pressure reduce speed in pneumatic system, mm of mercury/sec	2–5
Weight, kg	
AT-10, AT-11, AT-30, AT-31	0,45
AT-12, AT-13, AT-20, AT-21, AT-22	0,50
AT-41	0,80

## Safety information

---

Consult your doctor prior to using in pregnancy or if diagnosed with arrhythmia or arteriosclerosis. Please read this section carefully before using the unit.

- Always consult your doctor. Self-diagnosis of measurement results and selftreatment are dangerous.
  - People with severe blood flow problems, or blood disorders, should consult a doctor before using the unit, as cuff inflation can cause internal bleeding.
  - Do not leave the unit unattended with infants or persons who cannot express their consent.
  - Do not use the unit for any purpose other than measuring blood pressure.
  - Do not disassemble the unit or arm cuff.
  - Use only the approved arm cuff for this unit. Use of other arm cuffs may result in incorrect measurement results.
  - Make sure that the air tube is not wrapped around other parts of your body when taking measurements at night. This could result in injury when the air pressure in the air tube is increased.
  - Do not leave the cuff wrapped on the arm if taking measurements during the night. This could result in injury.
  - Do not inflate the arm cuff over 299 mmHg.
  - Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
  - Do not operate the unit in a moving vehicle (car, airplane).
  - Do not apply strong shocks and vibrations to or drop the unit and arm cuff.
  - Do not take measurements after bathing, drinking alcohol and smoking, exercising or eating.
  - Do not forcibly bend the arm cuff or bend the air tube excessively.
  - When removing the air tube, pull the air tube with holding near the air jack.
  - Do not expose the device to extreme temperatures, humidity, or direct sunlight.
- Never contact the cuff fabric with a sharp instrument, since this could cause damage.

## Direction for use

---

### Recommendations on blood pressure measuring

Measured blood pressure — blood pressure on vascular walls. It is generated by heart activity as a pump. At heartbeat pressure grows. Maximal pressure is called systolic (upper). During intervals between heartbeats blood pressure reduces and it is called diastolic (lower). To effect correct measuring one should know that blood pressure can vary abruptly even within short intervals. Blood pressure level depends on many factors. Generally it is lower in summer and higher in winter. Blood pressure changes depending on atmospheric pressure, physical activity, emotional lability, stresses and dietary habits. Medicines, digestions, alcohol, smoking are also important factors. Healthy people can have 30-50 mm of mercury measure difference of systolic pressure and up to 10mm of mercury measure difference of diastolic pressure. Blood pressure dependence on various factors is individual for each person. Thus, it is recommended to keep a diary of blood pressure indications.

Normal blood pressure for a healthy person should be 120/80mm of mercury. Anyway it is only a doctor who is competent to determine one's normal blood pressure and health condition.

#### 1. Environment

Blood pressure measure should be effected in quiet, calm and convenient conditions at comfort environmental temperature. Avoid external actions that possibly could increase measuring error or interfere measuring process and sit with your back straight nearby table.

The table height should provide blood pressure measuring in such a way that cuff midline is located at heart level, approximately at 4th intercostal space while in sitting position, or at axillary midline in reclining position.

Deviation in cuff midline position to heart level can result in incorrect blood pressure measure. Back support and hand placement on a supporting surface exclude blood pressure increase caused by muscles traction because of tension.

#### 2. Measuring preparation and rest time

Blood pressure should be measures 1–2 hours after meal. One should not smoke or take caffeine at least 1 hour before measuring. There should be no

tight clothes on a patient. Hand prepared for blood pressure measure should be bare. One should sit with his back leaned on a chair, legs rested. It is not recommended to talk during the measure, because it may influence the result. Repeated measure can be effected no sooner than in 5-minute rest.

### 3. Cuff size

Cuff air bladder length should be no less than 80% of arm circumference. Exact sizes of arm circumference are marked on a cuff. Cuff width should be about 40% of arm circumference (12–14 cm on average).

Generally blood pressure is measured on the right hand or on the hand with higher blood pressure level (when there are diseases resulting in significant blood pressure difference between right and left hand, as a rule the lower blood pressure level is indicated on the left hand). Usage of narrower or shorter cuff results in significantly wrong blood pressure level excess.

### 4. Cuff position

Detect brachial artery pulsation at arm midline level. Cuff bladder midline should be located exactly above the artery. Lower cuff edge should be located 2.5 cm above antecubital fossa (Fig. 2). Cuff tension: finger should pass between the cuff and arm surface.


Fig. 2

## 5. Determination of maximal air injection level

It is necessary for precise indication of systolic blood pressure and avoidance of auscultatory lapse:

- Detect radial artery pulsation, pulse rhythm and character. If there is distinct rhythm disturbance (fibrillation), systolic blood pressure level may vary from one contraction to another, and it is necessary to effect extra measuring to detect its exact level.
- Keeping feeling radial artery, one should quickly inflate the cuff pressure up to 60 mm of mercury, and then inflate 10 mm of mercury till pulsation disappears.
- Air should be let out of the cuff at 2 mm of mercury per sec. The registered blood pressure level is that when pulse appears again.
- Air should be totally let out of the cuff. To determine maximal air inflation level systolic blood pressure level should be increased by 30 mm of mercury.

## 6. Stethoscope position

The point of maximal brachial artery pulsation should be determined. Usually it is situated right above antecubital fossa inside arm.

Stethoscope diaphragm should completely fit to arm surface (Fig. 3). One should avoid too strong press of stethoscope because it may cause extra brachial artery compression. It is rec recommended to use low-frequency diaphragm.


Fig. 3

### **7. Cuff inflation and deflation**

Air inflation up to the maximal level should be effected quickly. Slow inflation may lead to venous outflow disturbance, strengthening of pain sensation and sound muting. Air should be let out of the cuff at 2 mm of mercury per second till Korotkoffs sounds appear. Then at 2 mm of mercury per second from sound to another. If the sound is too bad, air should be let out of the cuff quickly, the cuff position should be checked and procedure should be repeated. Slow deflation allows to determine systolic and diastolic pressure at the beginning of Korotkoffs sounds. Accuracy of blood pressure determination depends on decompression speed: the higher speed of decompression is the lower measuring accuracy is.

### **8. Systolic blood pressure**

Systolic blood pressure level is determined at Phase I of Korotkoffs sounds on the nearest scale mark (2 mm of mercury). If Phase I appears between two minimal scale marks, systolic pressure is considered to be the one corresponding to higher level. At significant rhythm disturbance additional blood pressure measuring is needed.

### **9. Diastolic blood pressure**

Diastolic blood pressure level is determined when the last distinct sound is heard. If Korotkoffs sounds continue till very low levels or zero, the registered blood pressure level corresponds to the beginning of Phase IV. Absence of Phase V can be observed among infants, expectant mothers and at states featured with high heart release level. In such cases diastolic blood pressure is determined as beginning of Phase IV of Korotkoffs sounds. If diastolic blood pressure is higher than 90 mm of mercury, measurement should be implemented within 40 mm of mercury. In other cases within 10–20 mm of mercury after the last sound disappears. Following this rule allows avoidance of determination of falsely high diastolic blood pressure at sounds recommencement after auscultator lapse.

### **10. Measurement results recording**

It is recommended to record which hand was used for measurement, cuff size and your position. The results are recorded as KI/KV. If Phase IV of Korotkoff's sound is determined the records are made as KI/KIV/KV. If there is no total disappearance of sounds, Phase V is considered to be equal to zero.

## 11. Repeated measurements

Repeated measurements should be implemented in 5 minutes after complete cuff deflation. Blood pressure level may vary from minute to minute. Average value of two or more measures made on one hand is usually more accurate than a single measure.

## 12. Blood pressure measure in other positions

It is recommended to measure blood pressure on both hands in lying and upright positions. Blood pressure changes are registered in 2–3 minutes after being in upright position. It is necessary to note which hand detects higher blood pressure level.

Blood pressure values of two hands may differ in more than 10 mm of mercury. The higher value more precisely corresponds to intra-arterial blood pressure.

### Special cases of blood pressure measuring

- "Auscultatory lapse" (temporary absence of sounds) may continue up to 40 mm of mercury. It is observed at high systolic blood pressure.
- "Infinite sounds" phenomenon is observed at high heart release: among infants, expectant mothers, at thyrotoxicosis, fever, aortic failure.
- Blood pressure measure among elderly persons. At elder age thickening and hardening of brachial artery wall is observed. It requires higher cuff pressure level (higher than intra-arterial pressure) to achieve thickened artery compression and it causes false overstating of blood pressure level (pseudo-hypertension phenomenon).
- Pulse reconnaissance on radial artery at cuff pressure level exceeding systolic blood pressure helps to determine blood pressure on forearm. If difference between systolic blood pressure values exceeds 15 mm of mercury actual blood pressure value can be determined by direct intra-arterial measuring. The patient should be informed of the problem occurred and corresponding record should be made in his case history to avoid possible mistakes in future.
- Extremely large arm circumference (adiposity, developed musculature), conical hand.
- Accurate blood pressure measuring can be impossible for patients having arm circumference over 41 cm or conical arm shape, because normal cuff position cannot be achieved. In such case it is recommended to use

corresponding cuff size and attempt to measure blood pressure level on upper arm.

### Phases of Korotkoff sounds

I phase	In this phase blood pressure is marked with steady sounds. Sound intensity increases gradually as cuff is released from air. The first of at least two sequential sounds is determined as systolic blood pressure.
II phase	Appearance of noise and rustling sound at further release of the cuff.
III phase	The period features crunch type increasing sound.
VI phase	This period corresponds to abrupt muting and appearance of soft blowing sound. This phase can be used for determination of diastolic blood pressure if the sounds can be heard till zero level pressure.
V phase	The period is characterized with elimination of the last sound and corresponds to diastolic blood pressure level.

### Environment

Storage temperature, °C	+10 – +35
Storage Humidity Range, %	20–85
Atmosphere pressure, hPa	700–1060

## Care and maintenance

### Disinfection

External surfaces of device and cuff can be treated with 3% peroxide solution mixed with 0.5% soap detergent solution of "Lotos" type or any other disinfectant in accordance with the applicable national standards.

If there is problems during measurement, check the following moments and take appropriate steps if necessary.

№	Defect	Possible cause and solution
1	Poor tone transmission, distortion or noise	<ol style="list-style-type: none"> <li>1. Check the integrity and conductivity of ear tips and binaural tubes. Be sure they are completely fit.</li> <li>2. Check the integrity and conductivity of the connecting tube. Make sure that junction between binaural tube and chest piece is sealed.</li> <li>3. Check the integrity and tightness of diaphragm to the chest piece.</li> </ol>
2	While inflating the cuff by the rubber balloon pressure is not increasing	<ol style="list-style-type: none"> <li>1. Make sure the drain valve is closed.</li> <li>2. Make sure that the cuff is correctly connected to a rubber bulb and manometer.</li> </ol>
3	Rate of outlet air can't be set at 2–3 mm of mercury/sec by adjusting of air outlet	<ol style="list-style-type: none"> <li>1. Disconnect the valve from the bulb in order to check whether there are obstacles inside the valve. Remove the obstacle and try again. If the valve is not working properly, replace it to avoid inaccurate measurement.</li> </ol>
4	While not using the device indicator is not at the level of $0 \pm 3$ mm of mercury	<ol style="list-style-type: none"> <li>1. Make sure that the valve is fully open before zero point set checking</li> <li>2. If the deviation from the zero point exceeds 3 mm of mercury, refer to the trade organization to recalibrate the monometer.</li> </ol>

---

**Attention!** In the event of other defects please contact your supplier.

---

### **Storage and transportation**

---

Transportation by suitable vehicles according to transportation requirements to prevent damage, deterioration or other adverse effects of medical devices.

### **Transport and storage conditions**

Temperature, °C	+10 – +35
Humidity Range, %	≤96
Atmosphere pressure, gPa	500–1060

### **Disposal information**

---

With reference to the disposal of the appliance, its accessories and its single use disposable parts, please see current waste disposal regulations.

Do not act directly on the appliance.

### **Warranty**

---

This appliance is warranted for a period of 12 months from the date of purchase as shown by the dealer's stamp and signature and by the purchaser's receipt for tax purposes which must be kept attached to this warranty.

■ Apexmed International B.V., Keizersgracht 62-64, 1015 CS, Amsterdam, the Netherlands. Ningbo Greatcare Trading Co., Ltd., Unit 93, Building 12, No. 818 Qiming Road, Yinzhou, 315105 Ningbo, Zhejiang, China (Китай).

**APEXMED**  
Healthcare essentials